

Polonaise In G Minor

Bach

Measures 1-4 of the Polonaise in G Minor. The piece is in 3/4 time and G minor. The first system shows measures 1 through 4. The right hand features a melodic line with slurs and fingerings (2, 5, 2, 2, 5, 1, 5, 4, 2, 5, 3, 4, 3, 1, 4, 5, 1, 2). The left hand provides a steady accompaniment with slurs and fingerings (4, 1, 4, 4, 1, 5, 5, 4, 2). The dynamic marking is *mf*.

Measures 5-8 of the Polonaise in G Minor. The second system shows measures 5 through 8. The right hand continues the melodic line with slurs and fingerings (2, 5, 2, 2, 5, 1, 5, 4, 2, 5, 3, 4, 3, 1, 4, 5, 1, 2). The left hand accompaniment has slurs and fingerings (1, 3, 2, 1, 4, 4, 1, 5, 1, 5, 5, 4, 2, 1, 3, 2, 5). The dynamic marking is *mp* for measures 5-7 and *f* for measure 8.

Measures 9-13 of the Polonaise in G Minor. The third system shows measures 9 through 13. The right hand features a melodic line with slurs and fingerings (2, 3, 1, 5, 4, 3, 4, 3, 1, 2, 1, 2, 2, 1). The left hand accompaniment has slurs and fingerings (3, 2, 5, 1, 2, 3, 5, 1, 3, 1). The dynamic marking is *f*.

Measures 14-16 of the Polonaise in G Minor. The fourth system shows measures 14 through 16. The right hand features a melodic line with slurs and fingerings (2, 1, 3, 1, 5, 4, 2, 1, 5, 4, 2, 1, 5, 1, 4, 3). The left hand accompaniment has slurs and fingerings (5, 2, 1, 2, 1, 5, 5, 1, 4, 3). The dynamic marking is *f*.

Measures 17-20 of the Polonaise in G Minor. The fifth system shows measures 17 through 20. The right hand features a melodic line with slurs and fingerings (1, 3, 1, 4, 1, 3, 1, 4, 1, 3, 1, 4, 4, 2). The left hand accompaniment has slurs and fingerings (2, 1, 2, 5, 2, 2, 1, 2, 5, 2, 3, 5, 5, 2, 2). The dynamic marking is *p*.

Measures 21-24 of the Polonaise in G Minor. The sixth system shows measures 21 through 24. The right hand continues the melodic line with slurs and fingerings (2, 5, 2, 2, 5, 1, 5, 4, 2, 5, 3, 4, 3, 1, 4, 5, 1, 2). The left hand accompaniment has slurs and fingerings (4, 1, 4, 4, 1, 5, 1, 5, 4, 2, 1, 5, 5, 4, 2, 1, 5). The dynamic marking is *mp*.